BD FACSLyric™

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Course Schedule

Day 1	Day 2	Day 3
• System basics • Maintenance • Universal setup • Compensation and settings	• Performance QC • Experiment setup • Worklist Acquisition and Analysis	• Data Management • Daily workflow • Settings practice • Fluidics Commands and Troubleshooting
Day 1 review	Day 2 review	• Final review

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What are the three main subsystems on the BD FACSLyric™?

Fluidics


Optics

Electronics

In BD FACSuite™ software, how many fluidic modes are available?


*For Research Use Only.

Where are the detector arrays?

What is the main purpose of the electronics subsystem?

Where can you perform quality control tasks in BD FACSuite?

Setup & QC

QC Reports

QC Tracking

Assay Setup Reports

Configurations

Current Configuration: 4-Blue 3-Red 3-Violet

Last Characterization QC (Normal Fluidics Mode): Passed Lot ID: 17884 1/26/2017 3:07 PM

Last Characterization QC (High Sensitivity Fluidics Mode): Passed Lot ID: 17884 1/26/2017 3:07 PM

Last Performance QC (Normal Fluidics Mode): Passed Lot ID: 17884 2/2/2017 10:35 AM

Last Performance QC (High Sensitivity Fluidics Mode): Passed Lot ID: 17884 2/2/2017 10:35 AM

Last LW Update: Passed Lot ID: 17884 1/30/2017 11:21 AM

Last LNW Update: Passed Lot ID: 17884 1/30/2017 11:21 AM

CS&T Bead Lot Transfer Summary

Configuration	Old CS&T...	New CS&T...	Completion...
4-Blue 3-Red 3-Violet			

Assay Setup Summary


Assay	Type	Status	CS&T Bead...	Completion...
6-Color Stimulation UD	User-Defined	OK	17884	2/2/2017 10:45 AM

Tube Settings Setup Summary


Tube Setting	Type	Status	CS&T Bead...	Completion...
6-Color Stimulation	User-Defined	OK	17884	2/2/2017 10:45 AM
Lyse No Wash	BD	OK	17884	2/2/2017 10:35 AM
Lyse Wash	BD	OK	17884	2/2/2017 10:35 AM

Which workspace allows you to import bead and reagent information?

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What is universal setup?

An automated system that establishes and maintains instrument MFI target values on a single instrument using...

BD FACSuite™ Software

BD™ CS&T Beads

BD™ FC Beads


Why is it important?

- Allows for reproducible and consistent results.
- Eliminates the need for daily compensation.
- Creates settings that are easily transported from instrument to instrument.

How it Works

Characterization QC is the initial setup task. Other tasks ensure all settings are updated.


CS&T beads are used during setup tasks to establish PMT voltages that will place cells at the same target everyday.

FC Beads are used to measure spillover values that are used in reference settings.

Universal Setup: Setup and QC Tasks

Characterization QC (CQC)

Beads

- Determines initial cytometer performance.

CS&T

Performance QC (PQC)

- Measures daily cytometer performance.

CS&T

Assay and Tube Settings Setup (ATSS)

- Updates user-defined settings based on daily performance.

CS&T


Update Reference Settings

- Re-measures spillover values.

CS&T
and
FC Beads

Universal Setup Workflow


Characterization QC

Characterization QC (CQC)

What you see...

What the software does...

- Performs laser alignment
- Defines initial %rCV, linearity, and sensitivity values
- Determines initial PMT voltages to place the bright beads at the following MFIs:
 - ~120,000 for scatter channels
 - ~100,000 for fluorescent channels

CQC Report

Cytometer Characterization QC Report 4-Blue 3-Red 3-Violet

Cytometer:	BD FACSLytic	User:	Admin User
Cytometer Name:	Julia	Institution:	None
Serial Number:	826591800039		
Fluidics Mode:	Normal		
Last Characterization & QC:	6/8/2017 2:51 PM		
Configuration Name: 4-Blue 3-Red 3-Violet		Last Modified: 4/5/2017 4:38 PM	

SUMMARY: PASSED

Target MFI and %rCV values

Linearity and sensitivity values

DETECTOR SETTINGS

Laser	Detector				PMT		Bright Bead		Mid Bead		Dim Bead		Linearity (±2%)			Resolution		
	Name	Minor	Filter	Position	Voltage	Slope of Gain	Median	%rCV	Median	%rCV	Median	%rCV	Min Channel	Max Channel	Sensitivity	Electron Noise RSD	Qr (x10 ³)	Rr
LASER: Blue (Wavelength = 488nm)																		
X	FSC	-	-	-	FSC	421.4	0.0040	120.533	0.6	120.685	0.6	10.912	0.6	N/A	N/A	568	N/A	N/A
X	SSC	10	488/15	E	442.0	8.3267	121.880	1.6	123.300	1.6	56.473	0.9	N/A	N/A	1.514	N/A	N/A	N/A
X	FTTC	507LP	527/32	D	616.8	8.6613	101.118	1.6	2.663	9.2	438	26.7	126	329.625	422	19.5	77.3	160
X	PE	560LP	565/42	C	549.9	8.5635	100.910	1.2	2.626	6.3	498	18.8	155	229.944	1.290	18.2	636.6	142
X	PerCP-Cy5.5	668LP	700/54	B	579.2	8.7237	99.647	2.1	2.262	9.5	404	22.6	84	220.892	512	20.2	36.2	18
X	PE-Cy7	752LP	783/56	A	625.0	8.9846	99.903	3.8	1.253	19.3	259	42.2	34	229.470	2.054	19.5	48.1	1
LASER: Red (Wavelength = 640nm)																		
X	APC	660/10	660/10	C	528.6	8.4813	101.205	1.7	2.805	8.0	494	19.6	120	226.756	568	19.3	52.5	23
X	APC-4700	705LP	730/40	B	507.8	8.5629	100.416	1.9	1.916	6.1	366	14.4	152	226.249	334	20.6	29.7	75
X	APC-Cy7	752LP	783/56	A	498.6	8.6896	101.237	2.4	1.642	6.1	337	13.4	91	224.520	323	19.9	41.3	138
LASER: Violet (Wavelength = 405nm)																		
X	V450	448/45	448/45	E	731.6	8.6550	100.828	2.6	6.076	8.6	747	37.1	153	228.177	153	21.0	72.3	1.253
X	V500-C	500LP	528/45	D	460.5	8.8811	101.807	2.3	2.641	7.2	341	21.7	95	228.462	120	21.8	39.0	977
X	RV505	006/36	006/36	C	693.9	8.4959	100.511	3.2	2.518	14.8	372	50.3	134	228.104	1.888	22.5	373.5	42

CQC Timeline

Performed by field service engineer (FSE)

Re-run by Admin user

- Every 6 months
- Using CS&T beads

Exercise

1. Import bead lots.

- a. Download bead lot file(s) from the website.
- b. Import file(s) into the library.

2. Run CQC.

- a. Prepare CS&T beads.
- b. Run CQC.
- c. Review results.

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
The screenshot shows a software window titled "Setup & QC Options". It contains the following elements:

- Task:** A dropdown menu currently showing "CS&T Bead Lot Transfer".
- Old CS&T Bead Lot ID:** A dropdown menu showing "46573 (Expires: 9/30/2016)".
- New CS&T Bead Lot ID:** A dropdown menu showing "86627 (Expires: 4/30/2016)".
- New Lot...:** A button located below the ID dropdowns.
- Start:** An orange button at the bottom left.
- Abort:** A grey button at the bottom right.

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


FITC Spillover


Compensation

The mathematical process by which the fluorescence values for the primary detector reflect only the fluorescence of the intended fluorochrome and the fluorescence from any other fluorochromes is removed.

Statistics	
Name	PE-A Median
FITC Stained Control: Unlabeled	55
FITC Stained Control: FITC	2,691


Statistics	
Name	PE-A Median
FITC Stained Control: Unlabeled	34
FITC Stained Control: FITC	35

Compensation Terminology


Undercompensated

Statistics	
Name	PE-A Median
Lyse Wash_FITC-FITC Negative	27
Lyse Wash_FITC-FITC Positive	707

FITC Negative	13
FITC Positive	3,119


Correctly compensated

Statistics	
Name	PE-A Median
Lyse Wash_FITC-FITC Negative	19
Lyse Wash_FITC-FITC Positive	19

FITC Negative	9
FITC Positive	9

Overcompensated

Statistics	
Name	PE-A Median
Lyse Wash_FITC-FITC Negative	17
Lyse Wash_FITC-FITC Positive	-193

FITC Negative	8
FITC Positive	-250


Tips for Accurate Compensation

- The fluorescence spectrum (% spillover) of the compensation control reagent should be identical to the reagent used in the experiment.
 - For example, do not use Alexa Fluor® 488 to compensate FITC.
- If a separate negative control tube is used, the cell/particle type must match that of the positive control tube.
 - For example, do not use lymphocytes for the negative tube and monocytes for the positive tube.
- Positive populations should be relatively bright and on scale.
- Collect enough events to get statistically accurate numbers.
 - In other words, do not use a rare marker as your compensation control.

Choosing a Compensation Control Type

Use cells if...

- Your cells are expressing a reporter protein, such as GFP.
- Your cells are stained with a fluorescent dye, such as PI or 7-AAD.
- The species of your antibody differs from available BD CompBead choices.


Use BD™ CompBeads if...

- You have low sample concentration.
- You are looking for antigens that are rare or variably expressed.


BD™ CompBead Particles


- Use the same antibodies as the experimental samples
- Conserve the experimental sample
- Create bright and uniform positive fluorescence peaks

Cells vs Beads

What is a tandem dye?




Tandem Dyes Are Special




Experiment: 3 cocktails, 3 different APC-Cy7 conjugates

CD8 CD3 CD4

Compensation control: CD3 APC-Cy7 used for all cocktails

Compensation control: lot-specific controls

Lot-Specific Compensation

Lot-specific compensation can be assigned when creating new reference settings or adding a fluorochrome to an existing setting.

Create Reference Settings

Select the controls and lot IDs that will be used to create Reference Settings.
Click Add or Delete to modify which tubes to run.
Acquisition will start when Next button is clicked.

CS&T Beads

CS&T Bead Lot ID: 17884 (Expires: 2/28/2020)

Control Tubes

Fluorochrome	Control Type	Label	Lot ID	Unstained
FITC	▼ FC Beads	▼ Generic	443331 (9/30/2018)	▼
PE	▼ FC Beads	▼ Generic	443332 (9/30/2018)	▼
PerCP-Cy5.5	▼ FC Beads	▼ Generic	667712 (12/31/2018)	▼
PE-Cy7	▼ CompBeads	▼ CD45RA	546873 (5/11/2018)	▼
APC	▼ FC Beads	▼ Generic	667711 (12/31/2018)	▼
APC-H7	▼ FC Beads	▼ Generic	6084967 (4/30/2018)	▼


How does compensation work in FACSuite?

- Setup and QC tasks re-calculate compensation based on daily cytometer performance.
- When PMT voltages (PMTVs) are adjusted, compensation values are re-calculated.
- Spillover values (SOVs) need to be re-measured every 60 days.

Classical SOVs and PMTV

Classical SOVs calculated in median fluorescence intensity (MFI) are PMTV dependent because MFI values are dependent on the PMTV.

New Terminology

Assigned BD Unit (ABD Unit)

- A BD standard for fluorescence that is PMTV independent.
- Each CS&T bead lot has defined ABD units for each channel.
- One ABD unit is defined as 1/40,000th of the fluorescence of a lymphocyte stained with anti-CD4 for a particular detector or fluorochrome.

Measured SOV (MSOV)

- Generated when reference settings are created or updated by running compensation controls.

Normalized SOV (NSOV)

- The actual spillover in ABD units.
- Determined by the overlap into optical filters.

Calculated SOV (CSOV)

- Calculated from existing measured NSOVs and tube settings.


NSOVs and PMTV

Normalized SOVs are calculated in ABD units.

$PMTV_{(FITC)} = 450V \quad 500V \quad 500V$

$PMTV_{(PE)} = 550V \quad 550V \quad 500V$


NSOVs calculated in ABD units are PMTV independent
because ABD units are PMTV independent.

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





BD-Defined

Start Here

Lyse/Wash and Lyse/No-Wash reference settings

User-Defined

Customize any or all of the following:

- Flow rate
- PMTV
- Threshold
- Area scaling
- Window extension
- Spillover values

Tube Settings

- A collection of values that place the positive population at the same position (MFI) whenever the tube settings are applied
- Do NOT contain spillover values
- Linked to lyse wash settings

Experiments using tube settings still have spillover values because they are calculated from a reference setting.

Tube Settings Contain...

- Tube target values (TTVs)
- Threshold value
- Flow rate
- Window extension
- Area scaling factor

NOTE: TTVs are only established for the parameters in the tube.

Tube Settings		
Created Date: 2/2/2017 10:45 AM		
Author: Admin User		
Short Description:		
Window Extension: 3		
Flow Rate: Medium		
Area Scaling Factor (ASF) Ratio		
FSC	1.00	
Violet	1.00	
Blue	1.00	
Red	1.00	

Parameter	Tube Target Value (TTV)	Threshold (Ov)
FSC	0.18	10,000
SSC	1.25	
FITC	0.22	
PE	0.95	
PerCP-Cy5.5	0.33	
PE-Cy7	0.73	
APC	0.24	
APC-H7	0.22	

Tube Target Values

- The ratio of median fluorescence intensity (MFI) to assigned BD unit (ABD) for a parameter
- Used to place CS&T beads at the target MFI during daily setup
- Instrument and bead lot independent
 - Transportable from one instrument to another.

$$\frac{\text{MFI}_{\text{CS\&T bright bead}}}{* \text{ABD}_{\text{PE}}} = \frac{2,000}{8,000} = 0.250 \text{ TTV}_{\text{PE}} \text{ for this assay}$$

*From the CS&T bead lot file

How it Works

1. Use your cells to adjust settings (flow rate, PMTV, threshold, etc.).
2. Create tube settings for your optimized tube using CS&T beads.
 - The MFI of the bright bead is recorded for each parameter in the tube.
 - TTVs are calculated by using the ratio of the MFI and the ABD value for each parameter.

Reference Settings

- A tube setting with associated measured SOVs
- Instrument dependent
- Updated during daily setup and QC tasks
- Spillover needs to be re-measured every 60 days

Tube Settings Spillover Values


Updated Date: 2/8/2017 1:20 PM Expiration Date: 4/9/2017 1:20 PM Remove Fluorochrome

Measurement: Area Height

Fluorochrome	FMTV	Spillover Values							
		RTC	PE	PerCP-Cy5.5	PE-Cy7	APC	APC-H7	BV421	BV510
		Generic	Generic	Generic	CD45RA 348273	Generic	Generic	Generic	Generic
RTC	642.5	100.00	1.18	0.03	0.17	0.01	0.06	0.02	1.87
PE	483.7	14.44	100.00	0.03	1.34	0.01	0.03	0.01	0.39
PerCP-Cy5.5	580.0	2.80	24.19	100.00	1.21	1.86	0.06	0.02	0.14
PE-Cy7	602.4	0.28	1.81	17.50	100.00	0.21	1.33	0.01	0.02
APC	551.0	0.00	0.01	1.28	0.03	100.00	2.05	0.00	0.00
APC-H7	553.9	0.00	0.00	6.64	9.12	13.96	100.00	0.00	0.00
BV421	536.5	0.17	0.03	0.08	0.03	0.03	0.19	100.00	1.73
BV510	487.7	14.04	1.51	0.30	0.17	0.10	0.54	0.75	100.00
FV557SV	518.9	0.61	13.18	0.05	0.14	0.06	0.06	0.90	4.95


Reference Settings Contain...

How it Works


1. Use your cells to adjust settings (flow rate, PMTV, threshold, etc.).
2. Create reference settings for your optimized tube using CS&T beads and appropriate compensation controls.
 - The MFI of the bright bead is recorded for each parameter in the tube.
 - TTVs are calculated by using the ratio of the MFI and the ABD value for each parameter.
 - Spillover values are measured using your compensation controls.

BD-Defined Reference Settings

Lyse/Wash (LW) and Lyse/No-Wash (LNW) Settings

- Created at installation by the field service engineer
- Contain tube settings for normal human lymphocytes
- Updated daily with Performance QC
- Spillover needs to be re-measured every 60 days

BD™ FC Beads

- Fluorescent beads that enable the software to determine spillover values.
- Each tube has a mixture of positive and negative beads.

Exercise

1. Prepare BD FC Beads.
2. Update the lyse/wash settings.

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

Update Reference Settings 對話方塊開啟。

Update Reference Settings

Select which tubes to update in the Run column. Ensure that the correct lot IDs are selected. Remove button will permanently remove compensation for the selected fluorochrome from Reference Settings. Acquisition will start when Next button is clicked.

Run	Fluorochrome	Control Type	Label	Lot ID	Remove Fluorochrome
☒	FITC	FC Beads	Generic	443331 (9/30/2018)	Part of FC Beads Kit Remove
☒	PE	FC Beads	Generic	443332 (9/30/2018)	Part of FC Beads Kit Remove
☒	PerCP	FC Beads	Generic	443333 (9/30/2018)	Part of FC Beads Kit Remove
☒	PerCP-Cy5.5	FC Beads	Generic	667712 (12/31/2018)	Part of FC Beads Kit Remove
☒	PE-Cy7	FC Beads	Generic	667713 (12/31/2018)	Part of FC Beads Kit Remove
☒	APC	FC Beads	Generic	667711 (12/31/2018)	Part of FC Beads Kit Remove
☒	APC-Cy7	FC Beads	Generic	443343 (9/14/2018)	Part of FC Beads Kit Remove
☒	APC-R700	FC Beads	Generic	5U3R700 (8/31/2017)	Part of FC Beads Kit Remove
☒	APC-H7	FC Beads	Generic	4160540 (7/31/2017)	Part of FC Beads Kit Remove
☒	V450	FC Beads	Generic	4160539 (6/30/2017)	Part of FC Beads Kit Remove
☒	V500-C	FC Beads	Generic	4160537 (7/31/2017)	Part of FC Beads Kit Remove
☒	BV605	FC Beads	Generic	5U3BV605 (9/30/201)	Part of FC Beads Kit Remove

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
Library		Tube Settings		
		Name	Created Date	Short Description
Assays		Lyse No Wash	1/25/2017	Lyse No Wash
Beads and Reagents		Lyse Wash	1/25/2017	Lyse Wash
Keywords				
Labels				
Tube Settings				
BD				
User-Defined				

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Run Performance QC (PQC) Daily to...

Measure and track cytometer performance settings.

Running PQC

What you see...

Setup Tasks

Normal Fluidics Mode

- ✓ Cytometer Initialization
- ✓ Laser Delays
- ✓ Beads Identification
- Determine Performance Values
- Laser Alignment Check
- Characterization Values Check
- Determine PMT Voltages

High Sensitivity Fluidics Mode

- Laser Delays
- Determine Performance Values

What the software does...

- Checks laser alignment
- Measures %rCV, linearity, resolution, and laser power
- Compares PMT voltages to values established during characterization QC
- Updates the spillover matrix and PMTVs for the LW and LNW settings

Viewing the PQC Reports

Cytometer Performance QC Report
4-Blue 3-Red 3-Violet

Cytometer:	BD FACStyrac-R	User:	Admin User
Cytometer Name:		Institution:	None
Serial Number:	11134		
Fluidics Mode:	Normal		
Last Characterization QC:	1/26/2017 3:07 PM		
Configuration Name:	4-Blue 3-Red 3-Violet	Last Modified:	2/8/2017 9:03 AM

SUMMARY: PASSED


%rCV, linearity, sensitivity

DETECTOR SETTINGS

Laser Info	Detector				PMTV		Bright Bead		Linearity (±2%)		Resolution			
	Name	Mirror	Filter	Position	Actual	Δ	Median	% rCV	Min Channel	Max Channel	Sensitivity		Qr (x10 ³)	Br
LASER: Blue (Wavelength = 488nm)														
X	FSC	-	-	FSC	429.5	-0.7	119.982	1.9	N/A	N/A	300	0	N/A	N/A
X	SSC	10	488/15	E	409.4	-0.9	119.938	3.5	N/A	N/A	300	0	N/A	N/A
X	FTTC	507LP	527/32	D	786.8	-1.8	99.906	2.8	41	230.023	300	0	41.8	83
X	PE	560LP	586/42	C	549.9	-1.2	99.966	2.3	36	230.017	300	0	377.3	244
X	PerCP-Cy5.5	665LP	700/54	B	638.2	-1.4	99.929	3.0	29	229.872	300	0	26.2	116
X	PE-Cy7	752LP	783/56	A	711.6	-1.4	100.058	5.4	40	229.930	300	0	28.3	9
LASER: Red (Wavelength = 640nm)														
X	APC	660/10	660/10	C	603.6	0.4	100.154	2.6	38	229.625	300	0	40.9	186
X	APC-R700	705LP	720/30	B	560.0	0.3	100.072	3.1	5	229.691	300	0	15.3	0
X	APC-Cy7	752LP	783/56	A	543.6	0.1	99.917	4.3	41	230.057	300	0	18.4	0
LASER: Violet (Wavelength = 405nm)														
X	V450	448/45	448/45	E	699.5	-1.5	99.963	3.5	26	229.752	300	0	99.9	1.839
X	V500-C	500LP	528/45	D	463.7	-0.9	100.144	3.1	28	229.820	300	0	58.7	897
X	BV505	606/36	606/36	C	708.6	-1.6	99.981	3.4	51	229.828	300	0	285.9	41

Tracking Performance Settings


Exercise

1. Prepare the CS&T research beads.
2. Verify the bead lot ID and select the **Performance QC** task.
3. Run performance QC.
4. View the performance QC reports (normal and high-sensitivity fluidic modes).
5. View Levey-Jennings charts in the **QC Tracking** tab.

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Workflow: Experiment Setup

Task	Location
Add any new fluorochromes to the configuration.	Setup & QC workspace
Add reagents and labels to the library.	Library
Add any new fluorochromes or tandems to the LW or LNW settings.	Setup & QC workspace
Set up the experiment workspace.	Experiment workspace
Optimize the settings.	Experiment workspace
Create tube settings.	Experiment workspace
Build an analysis report.	Experiment workspace
Create an assay.	Experiment workspace


6-Color Stimulation Experiment

- Sample type: human peripheral blood mononuclear cells (PBMCs)
- 2 tubes
 - Unstimulated
 - Stimulated: 6-hour incubation with SEB (staphylococcal enterotoxin B)
- Samples are fixed

Population Tree

6-Color PBMC Samples

- CD3 FITC
- CD69 PE
- CD45 PerCP-Cy5.5
- CD4 PE-Cy7
- CD154 APC
- CD8 APC-H7

Cytometer
4-Blue 3-Red 3-Violet

▼ Status

▲ PMT Voltages

Threshold Operation ☐ And ☒ Or Add Remove

Name	A	H	W	Voltage	Threshold
FSC	✓			223.9	10,000
SSC	✓			411.9	5,000
FITC	✓			643.6	5,000
PE	✓			468.5	5,000
PerCP-Cy5.5	✓			580.9	5,000
PE-Cy7	✓			603.3	5,000
APC	✓			530.9	5,000
APC-H7	✓			553.8	5,000

What makes these special?

Add Fluorochrome

To apply lot-specific compensation, first add the fluorochrome to the library and then to the LW settings.

New Reagent

Reagent Info

Product Name: CD4 PE-Cy7

Single Color: ☒

Fluorochrome: PE-Cy7

Label: CD4

1

Add labels and reagents to the library.

Setup & QC Options

Task: Add Fluorochromes

Reference Settings Name: Lyse Wash

CS&T bead lot ID: 17884 (Expires: 2/28/2020)

Start


2

Add fluorochromes to the Lyse Wash reference settings.

Tube Properties

Lot-specific compensation part two: in order to ensure that the correct spillover values are used, select a label and a lot ID.


Name	Label	Lot ID
FITC	CD3	
PE	CD64	
PerCP-Cy5.5	CD45	
PE-Cy7	CD4	67851 (11/10/2020)
APC	CD154	
APC-H7	CD8	46275 (10/18/2019)

What next?

Adjust settings in your experiment, as needed.

- Flow rate
- Threshold
- PMTV
- Window extension
- Area scaling factor


A blue star indicates that settings have been modified.


Threshold

Adjust the threshold to eliminate debris from the data.


No Threshold

Threshold on FSC 50,000


Threshold


Threshold - - - - -
Processed events ———
Events not processed ·····


Placing the Populations on Scale

PMTV too low

PMTV too high


Optimized PMTV

How do you save the changes?

Create Tube Settings

- CS&T beads are run to determine new target values
- Include tube target values (TTVs), thresholds, flow rate, window extension, and area scaling factors
- Can be recalculated daily to ensure comparable results from day to day and from system to system

Where are the settings saved?

Tube settings are categorized in the library as BD or User-Defined. The BD settings are the LW and LNW settings.


Library	Tube Settings					Delete	
	Name	Created Date	Author	Short Description	Reference Settings		
Assays	5-Color Stimulation	2/2/2017	Admin User			N	
Beads and Reagents	5-Color Stimulation.MI	2/28/2017	Admin User			N	
Keywords	IP Assay_HB	2/21/2017	Admin User			N	
Labels	Lab Practical_HB	2/9/2017	Admin User			N	
Tube Settings							
BD							
User-Defined							

If a tube setting has associated measured spillover values, (Reference Settings), then it is indicated in this column.

Where are the SOVs coming from?

Scenario: Applying a user-defined (UD) tube setting

Exercise

1. Add fluorochromes to the LW settings.
2. Set up the Experiment workspace.
3. Optimize settings.
4. Create tube settings.

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

The 'New Reagent' dialog box contains the following fields:

- Product Name:** CD4 PE-Cy7
- Single Color:** ☒
- Fluorochrome:** PE-Cy7 (dropdown menu)
- Label:** CD4 (dropdown menu)

The 'Add New Lot' dialog box contains the following fields:


- Lot ID:** 67851
- Expiration Date:** 11/10/2020 (calendar icon)
- Buttons:** OK, Cancel

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Add Fluorochromes


Click Add or Delete to modify which tubes to run.
Acquisition will start when Next button is clicked.

Fluorochrome	Control Type	Label	Lot ID	Unstained
PE-Cy7	CompBeads	CD4	67851 (11/10/2020)	
APC-H7	CompBeads	CD8	46275 (10/18/2019)	

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Lyse Wash Print

Tube Settings **Spillover Values**

Updated Date: 3/27/2017 3:25 PM Expiration Date: 3/31/2017 11:21 AM Remove Fluorochrome

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1. Build a report.
2. Create an assay.

Building a Report

- Create and format plots.
- Adjust gate properties.
- Analyze data from multiple tubes.
- Add statistics to the report.

Worksheets and Reports

Worksheets

- Contain plots, gates and statistics
- Can include headers and footers

Reports


- Contain plots, gates and statistics
- Can include headers and footers
- Can be automatically exported
- Are associated with assays
- Can display an E-Signature


Plots and Gates


Dot Plot


Density Plot


Contour Plot

Histogram

Gate Properties




- Modify the name and color of a gate
- Display outlines, statistics or labels


Strict Parameter Matching

When strict parameter matching is selected, the long parameter names must match in order for the gate to be displayed.

Unique Gates



The unique gate feature allows gates to be tube specific. This allows changes to be made to a gate on one tube, but not affect the same gate on another tube.

Gate Hierarchy

The hierarchy shows the name, color and relationships with other gates.

Quadrants are considered statistical gates. They can be displayed by selecting the checkbox.


Statistics

- Statistics can be displayed for multiple tubes or can dynamically change based on the run pointer.
- Expressions and keywords can be added to statistics views.
- Among other attributes, the statistics properties can allow for the automatic export of statistics.

Create an Assay


Assays

- Contain a set number of tubes that share the same tube settings, worksheets, and reports.
- Are run by worklists using the Loader.
- Have elements that can be exported automatically.

Edit Assay Properties

Assay properties control which results are sent to the laboratory information system (LIS), which reports are printed or exported, E-Signature properties, and automatic approval options, among other options.

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
Create Assay

Enter an assay name with a maximum of 27 characters. UD (for user defined) will be added to the end of the name when it is saved.

Name: 6-Color Stimulation

Description: 2 tubes: Unstimulated and Stimulated

Report to Display After Entry Run: Analysis

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

Select Select All Remove Remove All Move Up Move Down

Selected

- All Events, lymphocytes, CD45 +, CD4+, Q2, Q2_UR, % Parent (unstimulated)
- All Events, Events (unstimulated)
- All Events, lymphocytes, CD45 +, CD4+, Q2, Q2_UR, % Parent (stimulated)
- All Events, % Parent (unstimulated)
- All Events, % Grandparent (unstimulated)
- All Events, % Total (unstimulated)
- All Events, lymphocytes, CD45 +, CD8+ T cells, CD8+ T cells_CD3- CD8+, Events (unstimulated)
- All Events, lymphocytes, CD45 +, CD8+ T cells, CD8+ T cells_CD3- CD8+, % Parent (unstimulated)
- All Events, lymphocytes, CD45 +, CD8+ T cells, CD8+ T cells_CD3- CD8+, % Grandparent (unstimulated)

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Select Select All Remove Remove All

Selected

- All Events, Lymphocytes, CD45+, CD4 T-cells, Q2, Q2_UR, % Parent (Unstimulated)
- All Events, Lymphocytes, CD45+, CD4 T-cells, Q2, Q2_UR, % Parent (Stimulated)

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Before You Begin

For User-Defined assays

1. Create tube settings.
2. Create an assay.
 - Can contain multiple tubes and tube settings
 - Create a report
 - Set assay properties

Worklist Preferences


Acquisition

- Acquisition Delay Timer
- Report Delay Timer

Export


- Entry Run Package
- FCS
- Result

Print

Building a Worklist

1. Select a loading option.

2. Select tasks to be run.

3. Enter sample IDs.

Worklist Entries					
	Sample ID	Task	Status	Location	Sample Carrier
# 1	Donor 1	6-color Stimulation UD	Ready	A1-A2	001
1.1		Unstimulated	Ready To Acquire	A1	001
1.2		Stimulated	Ready To Acquire	A2	001

4. Rename the worklist.


Acquiring Data

1. Edit carrier properties (optional).
2. Load the carrier.
3. Click **Run All**.
4. Verify gates and PMT voltages.
 - Use the stop timer, if necessary.


Analyzing Data


1. Verify gates and adjust them, if necessary.
2. Print the report.

Putting the Pieces Together

Exercise

1. Acquire data.


- Set worklist preferences.
- Build the worklist.
- Run samples.

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Acquisition Delay Timer

Preview for 20 seconds

☒ Use timer for the first tube in each entry (audible alarm will sound)

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Auto Export Options

☒ Export entry run package after acquisition

Location

C:\BD Export\ERP\Worklists

Browse ...

☒ Create dated subfolder (yyyyMMdd)

Naming Format

Worklist Name

Sample ID

None

None

Delimiter

Underscore

☒ Autonumber starting with

1

Example based on the selected choices

Worklist Name_Sample ID_001.AssayEntryRun

FCS

Auto Export Options

☒ Export FCS after acquisition

Location

C:\BD Export\FCS\Worklists

Browse ...

☒ Create dated subfolder (yyyyMMdd)

Naming Format

Worklist Name

Sample ID

Tube Name

None

Delimiter

Underscore

Example based on the selected choices

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
Worklist Entries				
Approval Audit Trail E-Signature More...				
	Sample ID	Task	Status	Location
▶ 1	Customer 1	6-color Stimul.	Ready	-
▶ 2		Perform Daily	Ready	-

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Worklist Entries							
Approval Audit Trail E-Signature More...							
	Sample ID	Task	Status	Location	Acquisition Time	Tube ID	Case Number
▶ 1	Customer 1	6-color Stimul.	Ready	-			DP1109
▶ 2		Perform Daily	Ready	-			

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The screenshot shows the 'Auto Export Options' dialog box. Under 'Export entry run package after acquisition', the 'Location' is set to 'C:\BD Export\ERP\Worklists' and 'Create dated subfolder (yyyyMMdd)' is checked. In the 'Naming Format' section, the dropdowns are set to 'Worklist Name', 'Sample ID', 'None', and 'None'. The 'Delimiter' is 'Underscore' and 'Autonumber starting with' is '1'. An example filename is shown as 'Worklist Name_Sample ID_001.AssayEntryRun'.

FCS

The screenshot shows the 'Auto Export Options' dialog box. Under 'Export FCS after acquisition', the 'Location' is set to 'C:\BD Export\FCS\Worklists' and 'Create dated subfolder (yyyyMMdd)' is checked. In the 'Naming Format' section, the dropdowns are set to 'Worklist Name', 'Sample ID', 'Tube Name', and 'None'. The 'Delimiter' is 'Underscore'. An example filename is shown as 'Worklist Name_Sample ID_Tube Name_yyyyymmdd_hhmmss.fcs'.

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The screenshot shows the 'Auto Export Options' dialog box. Under 'Automatically export results', the 'Location' is set to 'C:\BD Export\Results\Worklists' and 'Create dated subfolder (yyyyMMdd)' is checked.

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
Worklist Entries					
Approval Audit Trail E-Signature More...					
	Sample ID	Task	Status	Location	Carrier ID
▶ 1	Customer 1	6-color Stimulation	Ready	A1-A2	001
▶ 2		Perform Daily Clean	Ready	A3-A4	001

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Worklist Entries								
Approval Audit Trail E-Signature More...								
	Sample ID	Task	Status	Location	Carrier ID	Acquisition Time	Tube ID	Case Number
▶ 1	Customer 1	6-color Stimul.	Ready	A1-A2	001			DP1109
▶ 2		Perform Daily	Ready	A3-A4	001			

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Fluidics Commands

SIT Flush

Backflushes sheath fluid through the SIT and removes bubbles from the flow cell.

Clean Cuvette

Runs cleaning solution through the flow cell.



Drain and Fill Flow Cell

Removes persistent bubbles from the flow cell.

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

Assay Setup Report
9-Color IP_HP UD

Cytometer Model	BD FACSLyric-R	Bead Product	BD CS&T Beads	User	Admin User
Cytometer Name		Lot ID	17884	Institution	None
Serial Number	L1234	Expiration Date	2/28/2020		
Configuration Name	4-Blue 3-Red 3-Violet				


SUMMARY: PASSED

WARNINGS

Spillover values for one or more fluorochromes are not established

Associated reference settings expired

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
Worklist Entries					Approval	Audit Trail	E-Signature	More...
	Sample ID	Task	Status	Location				
1	Customer 1	6-color Stimulation_MP UD	Ready For Approval	A1-A2				
1.1		Unstimulated	Complete	A1				
1.2		Stimulated	Complete	A2				
2	Customer 2	6-color Stimulation_MP UD	Ready	A3-A4				
2.1		Unstimulated	Ready To Acquire	A3				
2.2		Stimulated	Ready To Acquire	A4				
3		Perform Daily Cleaning	Ready	A5-A6				

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
Daily Workflow: Setup and QC


- Run Performance QC.
- Run Assay and Tube Settings Setup, if needed.

Setting up the System


Assay and tube settings setup updates user-defined tube settings and assays to reflect the current cytometer performance.

Recalculating Tube Settings


Assay and Tube Settings Setup

Voltages are adjusted to place the bright CS&T beads at the saved tube target values.


Assay and Tube Settings Setup

Select all assays or tube settings you would like to run that day.

Running Assay and Tube Settings Setup

What you see...

What the software does...

Uses CS&T beads to determine PMT voltages needed to reach user-defined settings.

Assay Setup Report

Assay Setup Report 6-Color Stimulation UD

Cytometer Model	BD FACSymyx-R	Bead Product	BD CS&T Beads	User	Admin User
Cytometer Name		Lot ID	17884	Institution	None
Serial Number	L1234	Expiration Date	2/28/2020		
Configuration Name	4-Blue 3-Red 3-Violet				

SUMMARY: **PASSED**

PMTV and Spillover Values - Stimulated (6-Color Stimulation)

Parameter	PMTV	Bright Bead Median		Spillover Values					
		Target	Actual	FTTC Generic	PE Generic	PerCP-Cy5.5 Generic	PE-Cy7 Generic	APC Generic	APC-H7 Generic
FSC	223.9	17,986	17,987	N/A	N/A	N/A	N/A	N/A	N/A
SSC	431.9	125,232	125,104	N/A	N/A	N/A	N/A	N/A	N/A
FTTC	643.6	16,983	16,955	100.00	1.13	0.03	0.17	0.01	0.06
PE	468.5	24,267	24,225	13.31	100.00	0.03	1.47	0.01	0.03
PerCP-Cy5.5	198.9	43,242	43,190	3.12	23.18	100.00	1.13	2.00	0.09
PE-Cy7	605.1	23,251	23,210	0.29	1.77	16.60	100.00	0.29	1.44
APC	530.9	32,199	32,631	0.00	0.01	0.89	0.02	100.00	0.95
APC-H7	553.8	117,659	117,799	0.00	0.00	5.26	7.49	19.10	100.00

Exercise

1. Run performance QC.
2. Verify the CS&T bead lot ID and select the **Assay & Tube Settings Setup** task.
3. Select the assay(s) and/or tube settings you are running.
4. Run assay and tube settings setup.
5. View the **Assay Setup Report**.

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

Are all of the fluorochromes in
the configuration?

Where would you check?

*BD Horizon™ Fixable Viability Stain

Experiment Setup: SOVs

1. Do all of these fluorochromes exist in the LW spillover matrix?

Yes

Great! Proceed to question 2.

No

Add fluorochromes to the LW settings, then proceed to question 2.

2. Is the compensation adequate for your assay?

Yes

Continue using the LW settings.

No


Create reference settings.


Checking Compensation

Potential Compensation Errors

Correctly Compensated

Exercise

1. Choose a lab practice option.
2. Add a fluorochrome to the configuration, and labels and reagents to the library.
3. Continue with your option.
4. Analyze data and create an assay.
5. More settings practice.

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
CD117	BD
FMC-7	BD
User-Defined Label_001	FACSuite Student

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
Reagent Info	
Product Name	CD45RA PE-Cy7
Single Color	☒
Fluorochrome	PE-Cy7
Label	CD45RA

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Enter lot information for CD45RA PE-Cy7.

Lot ID	58612
Expiration Date	3/1/2017

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

The 'Setup & QC Options' dialog box shows the following configuration:

- Task:** Add Fluorochromes
- Reference Settings Name:** Lyse Wash
- CS&T bead lot ID:** 17884 (Expires: 2/28/2020)
- Start** button is highlighted in orange.

The 'Add Fluorochromes' dialog box contains the following information:

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Fluorochrome	Control Type	Label	Lot ID	Unstained
PE-Cy7	CompBeads	CD45RA	546873 (5/11/2018)	

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9-Color IP_HP UD Edit

General Export Results Reports Send to LIS Assign Lot ID

Assign default Lot ID to use lot-specific compensation for the tube.

	FITC	PE	PerCP-Cy5.5	PE-Cy7	APC	APC-H7	BV421	BV510	FVS575
	Assign Lot	Assign Lot	Assign Lot	Assign Lot	Assign Lot	Assign Lot	Assign Lot	Assign Lot	Assign Lot
Unstained	CD3 Generic	CD16+56 Generic	CD45 Generic	CD45RA Lot 546873	CD19 Generic	CD8 Lot 46275	CD45RO Generic	CD4 Generic	Viability Generic

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
The screenshot shows the 'Header/Footer' dialog box with two main sections: 'Header' and 'Footer'. Each section has a 'Show Border' checkbox checked. The 'Header' section contains two text boxes: the top one contains '<Worksheet/Report Name Label>' and '<Worksheet/Report Name>', and the bottom one is empty. The 'Footer' section contains two text boxes: the left one contains '<Experiment/Assay Name Label>' and '<Experiment/Assay Name>', and the right one contains 'Page 1 of N'. Red boxes highlight the text boxes in both sections and the 'Show Border' checkboxes.

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
Standardization


Standardization from Day to Day

Goal: To be able to compare results across multiple days and reduce user-to-user variation.

Standardization Across Cytometers

Goal: To be able to compare results across sites and cytometers.

Standardized Setup with BD FACSuite Software

One Cytometer—Day to Day

- Run performance QC daily.
- Create tube settings for all experiments and assays.
- Run Assay and Tube Settings Setup.

Multiple Cytometers

- Run performance QC on each cytometer.
- Import assays, experiments, or tube settings.
- Run Assay and Tube Settings Setup.

Transferring Assays

1. Run performance QC.
2. Create an assay with tube settings.
3. Export the assay from the library.
4. Import the assay into the library.
5. Run performance QC and assay and tube settings setup.
6. Run the assay.

Tips

- Import any reagents from the library as well.
- Verify that you have added all new fluorochromes to the configuration and reference settings.
- You can also import experiments and tube settings on other BD FACSLyric systems.

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© 2018 年，Becton, Dickinson and Company。保留所有權利。在未事先取得 BD Biosciences 書面許可下，不得複製、傳送、轉載或將本刊物的任一部分儲存於可叫出的系統或是以任何形式或任何方式翻譯成任何語言或電腦語言：電子、機械、電磁、光學、化學、手冊等。