

The New Milestone of Confocal Microscopy

- New full range confocal platform -







Drosophila melanogaster, eye section Red: F-Actin, Cy3 Blue: Nuclei, DAPI Green: pigmented cells, GFP Courtesy of Anne Galy, IGBMC, Strasbourg-Illkirch, France

Confocal Image -- Fluorescence

Confocal Microscopes can optically remove all information that is from outside the depth of focus.



Optical Section

The consequence is a sharp optical section.













Objective	Immersion	coverglass
HC PL APO 10x/0.4 CS	Air	0.17mm
HC PL APO 20x/0.7	Air	0.17mm
HC PL APO 40x/1.30 Oil CS2	Oil	0.17mm
HC PL APO 63x/1.40 Oil CS2	Oil	0.17mm
HC PL APO 100x/1.4 Oil CS2	Oil	0.17mm





- CS: Confocal Scan, the highest quality objectives
- New set of objective lenses for UV/VIS
- Optimized for multi-color overlay between 405 nm and VIS lasers
- Excellent planarity of field











Leica TCS SP8

Looking forward to your discoveries

- Sensitivity by Design
- Faster for new Biology
- New technology for image
- Ready to grow



Leica TCS SP8 Sensitivity by design



New scan optics

- New coatings for improved transmission
- Dedicated scan optics maximizing photon efficiency

 optimal transmission for each application
- Arranged in "4f" design for even field illumination
- "My Confocal" dedicated intermediate optics





New scan optics Coatings of scan optics compared 100% T[%] T[%] after passage of 8 per opt. 96% surface opt. surfaces 92% $(96\%)^8 = 72\%$ 96% No coating 88% **Transmissio** $(99.3\%)^8 = 95\%$ 99.3% (400-1300nm) 84% 80% 76% Reduced reflection by 60%72% VISIR HIVIS No coating Transmission after passage of 8 optical surfaces 72% 95% 98%





- Highly transparent → Increase sample lifetime, cell viability → more and better experiments
- Flexible → links in with WLL, Ar laser
 → more diverse experiments
- Ease of use → Save training time; better turnover (low failure rate)
- **Supports new dyes** → Future-proof
- Cool technology \rightarrow Prestige



AOBS Principle







Confocal microscopy: sophisticated

Leica Innovation Goals: Introduction of tuneable optical elements



Leading in Multispectral Imaging: Spectral Imaging Detector SP





- All dyes: Freely tuneable emission bands
- Low sample photodamage: high efficiency
- Up to 5 true confocal channels simultaneously
- Intuitive operation



Spectral scans

No distortions in spectra caused by filter or splitter transmission spectra





- Laser Kit White light laser
- 470-670nm
- WLL now supports
 - Lambda²Excitation-Emission scans
 - Lightgate (reflection removal)
- Add optional UV lasers
- LightGate: Non-optical reflection removal
- Pulsepicker option for improved FLIM
- Lifetime Saver











共軛焦顯微鏡中白光雷射的應用

✓ 只有搭配在AOBS系統上,才能真正稱之為白光雷射系統。







Leica HyD



Reduced light dosage: Increased cell viability



Improved Signal-to-Noise



Improved contrast of HyD vs. PMT





Faster for new biology



Single Point Scan Scanner





- Better sampling at low magnification
 - Large FOV
 - \circ Field number = 22
 - New Zoom range 0.75x .. 48x
- Save time on post-processing (mosaics)
- Homogenous illumination -> better data
- Now with 3600 Hz scan frequency (bidir)
- Equiv. to 7 fps @ 512 x 512







New technology for image

 $\langle \lambda 2 \rangle$ -Mapping





λ^2 –Mapping: Discover hidden information





Benefits of λ^2 -Maps

- Full spectral analysis of images
- Understand samples with very complex fluorescence
- Depict multifaceted features at one look
- Complete spectrum in each pixel
- Extract excitation and emission spectra



ADVANCED / USER CASE / APPLICATIONS

APPLICATION: FIND THE BEST FLUORESCENT LABEL

λ²–Mapping: **Find the best fluorescent label!**



Overview image



 λ^2 plots: Excitation: 470 to 670 nm, Detection: 500 to 750 nm



Strong autofluorescence of primary mouse hepatocytes

Sample: Courtesy of René Meyer, Klingmüller Group, Systems Biology of Signal Transduction, DKFZ, Heidelberg, Germany



ADVANCED / USER CASE / APPLICATIONS

APPLICATION: OPTIMIZATION OF DATA QUALITY

$\lambda^2-Mapping:$ Experiment optimization for best data quality





Sample: Mixture of fixed cells expressing four different fluorescent proteins. Excitation emission peaks are given in brackets. All cells show a small autofluorescence peak at (512, 533) nm.

Courtesy of Kees Jalink, Department of Cell biology, The Netherlands Cancer Institute Amsterdam, The Netherlands





45678 ON



80 MHz



<u>Sample:</u> Mixture of fixed cells expressing four different fluorescent proteins. XI Port: Mirror Wirror Action Standard Construction Standard Const

Courtesy of Kees Jalink, Department of Cell biology, The Netherlands Cancer Institute Amsterdam, The Netherlands



Major Instruments Co., Ltd.

New technology for image

Light Gate imaging



- The WLL is a pulsed laser
- Use time information to
 - \Rightarrow separate fluorescence decay from reflection
 - ⇒ discriminate between wanted and non-wanted fluorescence

Fluorescence lifetime Average time that molecules stay in their excited state



TCS SP8X LightGate Filter free removal of non-wanted background



- Detector data reading is switched off during WLL pulse.
 ⇒ Only fluorescence signal is used for image formation, reflected light is excluded.
- LightGate is adjustable to different fluorescence decay times.
 ⇒ Removal of non-wanted fluorescence is light possible.



Removal of reflected and scattered light by LightGate



Sequential scan with and without LightGate. Detection directly at excitation wavelength.



Removal of reflected and scattered light by LightGate



Excitation: 470 nm Excitation: 510 nm Detection: 465-590nmDetection: 495-540 nm C

Overlay



Sample: fixed HeLa cells, tubulin stained with BD Horizon V-500, nucleus stained with Chromeo 505, xy scan







New Software LAS X





Multidimensional Confocal Imaging ху



Drosophila leg, FITC, non-confocal

confocal

ху













xyzt







Fish embryo (Medaka) Flow of red blood cells and migration of macrophages. Macrophage: VFP and RFP (1st & 2nd ch) Red blood cells: TLD

Courtesy of Clemens Grabher and Jochen Wittbrodt (EMBL), Heidelberg, Germany

Multidimensional Confocal Imaging

xyzt



Drosophila mitosis pre-cellularized embryo GFP-tagged histone,

400 frames, 20 sec between each stack, each frame presents a stereo image

Courtesy: Prof. Sullivan (University of California, Santa Cruz); Robert Saint (Adelaide, Australia)

xyzt







共軛焦顯微鏡的發展趨勢:

Confocal system is not only for confocal images, but a platform for new application development...





