

High-Throughput Gene Expression Research

王簾讀

Kevin Wang


源資國際生物科技股份有限公司

Tri-I Biotech, Inc.

www.tri-ibiotech.com.tw



4 Stage Segmentation for Gene Expression Research



Customers in gene expression research fit into one or more of these segments.

Each segment has its own needs.

SuperArray offers products to meet the customer's need in each segment.




4 Stage Segmentation for a typical Gene Expression Research project.

- **1 - Discovery – High density genome wide**
 - Low sample #, High gene content
 - Affymetrix Services – Agilent?
- **2 – Screening – Low Density Pathway focused**
 - High sample number, lower gene number
 - SuperArray Oligo, RT PCR, cDNA Arrays
- **3 – Confirmation – RT PCR**
 - Very low gene number, high value genes
 - Superarray PCR products
- **4- Confirmation – Translation level**
 - Functional Biology
 - CASE Kits
 - siRNA




RT²Profiler™ PCR Array: Data Analysis & Tech Support




SuperArray Bioscience Corporation



HOW THE PCR ARRAY WORKS:



HOW THE PCR ARRAY SYSTEM 2.0 WORKS:



Post-Run: Baseline and Threshold

- ❖ **Manual Setting**
- ❖ **Baseline**
 - ❖ Use Linear View
 - ❖ Set to Cycle #2 up to 1 or 2 cycle values before earliest amplification
- ❖ **Threshold Value**
 - ❖ Use Log View
 - ❖ Place in
 - 1) Linear phase of amplification curve
 - 2) Above background signal, but within lower third of curve
- ❖ **Human 18S rRNA Check**
 - ❖ Omit if C_t value < 10
- ❖ **Export C_t values to blank spread sheet (Excel).**
- ❖ **Baseline/Threshold Must Be Same Between Runs**



Baseline and Threshold: Definitions

Baseline Value

During PCR, changing reaction conditions and environment can influence fluorescence. The background signal is most evident during the initial cycles of PCR prior to significant accumulation of the target amplicon. During these early PCR cycles, background signal in all wells is used to determine the “baseline fluorescence” across the entire reaction plate. The goal of data analysis is to determine when target amplification is sufficiently above the background signal, facilitating more accurate measurement of fluorescence.

Threshold

The threshold is the numerical value assigned for each run, statistically significant point above the calculated baseline.

Ct Value

The Threshold Cycle (Ct) reflects the cycle number at which the fluorescence generated within a reaction crosses the threshold.

From AppliedBiosystems Technical Note



Data File: APH_024_UniversailcDNA.opd

SYBR-490
Select a Reporter

Select analysis mode: PCR Base Line Subtracted Curve Fit

Select Wells
Reports

Threshold Cycle Calculation


Baseline Cycles
2 through 15
 Auto Calculated
 User Defined

Threshold Position
22.6
 Auto Calculated
 User Defined

Recalculate Threshold Cycles

Normal View

Save for X-axis Allelic Analysis
Save for Y-axis Allelic Analysis



Save data from current dye layer for allelic analysis (Y-axis).

Well	Cycle Ct
A1	19.8
A2	26.5
A3	23.8
A4	25.3
A5	26.4
A6	21.7
A7	31.0
A8	27.0
A9	23.5
A10	21.4
A11	23.8
A12	23.9
B1	26.7
B2	25.1
B3	27.1
B4	27.0
B5	25.6

Library

Workshop

Run-Time Central

Data Analysis

Running: Instrument not in Remote

9/15/2006 11:28 AM

Fold Change Calculation

Due to the inverse proportional relationship between threshold cycle and the original gene expression level, and the doubling of the amount of product of every cycle, the original expression level (L) for each gene of interest is expressed as:

$$L = 2^{-C_t}$$

To normalize the expression level of a gene of interest (GOI) to a housekeeping gene (HKG), the expression levels of the two genes are divided:

$$\frac{2^{-C_t(\text{GOI})}}{2^{-C_t(\text{HKG})}} = 2^{-[C_t(\text{GOI}) - C_t(\text{HKG})]} = 2^{-\Delta C_t}$$

To determine fold change in gene expression, the normalized expression of the GOI in the experimental sample is divided by the normalized expression of the same GOI in the control sample:

$$\frac{2^{-\Delta C_t \text{ expt}}}{2^{-\Delta C_t \text{ control}}} = 2^{-\Delta \Delta C_t} \quad \text{Where } \Delta \Delta C_t \text{ is equal to } \Delta C_t \text{ expt} - \Delta C_t \text{ control}$$

The complete calculation is as follows:

$$\frac{\frac{2^{-C_t(\text{GOI}) \text{ expt}}}{2^{-C_t(\text{HKG}) \text{ expt}}}}{\frac{2^{-C_t(\text{GOI}) \text{ control}}}{2^{-C_t(\text{HKG}) \text{ control}}}} = \frac{2^{-[C_t(\text{GOI}) - C_t(\text{HKG})] \text{ expt}}}{2^{-[C_t(\text{GOI}) - C_t(\text{HKG})] \text{ control}}} = \frac{2^{-\Delta C_t \text{ expt}}}{2^{-\Delta C_t \text{ control}}} = 2^{-\Delta \Delta C_t}$$



PCR Array System 2.0



Technical Sales Training



RT²Profiler PCR Array

❖ Pathway Focused

Profile the expression of a panel of genes relevant to a pathway or disease state.

❖ Simple and Accurate

Simple real-time PCR method provides high sensitivity and wider dynamic range.

❖ Designed for Routine Use

Brings expression profiling to almost any lab with a real-time PCR instrument.

❖ Continued successful marketing message

❖ Today: 2.0 = higher-level customer satisfaction



ISSUES AND PROBLEMS

❖ RNA Quality Control

- ❖ Purity and Integrity, Genomic DNA Contamination

❖ Sensitivity

- ❖ Lower (ng) RNA amounts
- ❖ Lower C_t values, higher positive calls

❖ Convenience

- ❖ Too many cocktails and dilutions, esp. for 384-well

❖ PCR Efficiency and Consistency

- ❖ Internal PCR control, Replicates of Controls



SOLUTION = PCR ARRAY SYSTEM 2.0

❖ RT²Profiler PCR Array

❖ With NEW Control Features & Simpler Protocol

Internal & External Control, gDNA Detection, Only One Cocktail

❖ RT² Real-Time SYBR Green PCR Master Mixes

❖ NEW: RT² PCR Array First Stand Synthesis Kit


Improves and Controls for RT Efficiency

❖ NEW: RT² PCR Array RNA QC Plate


Tests for genomic DNA and inhibitors before running PCR Array



HOW THE PCR ARRAY WORKS:



HOW THE PCR ARRAY SYSTEM 2.0 WORKS:




NEW PCR ARRAY AVAILABILITY

- ❖ **40 Pathways available for Human & Mouse**
 - ❖ cDNA GEArray and MultiGene-12 conversions
- ❖ **12 Pathways for Rat**
- ❖ **96-well (A,B,C,D) and 384-well (E) formats**
- ❖ **Customization**



New Format: 384-well PCR Array (Type E)

- Specific for customers using ABI 7900HT with 384 well block
- Sold in sets of 4 plates



RT² PCR Array RNA QC Plate LAYOUT

HK1	HK1	HK1	HK1	HK1	HK1	HK1	HK1	HK1	HK1	HK1	HK1
HK2	HK2	HK2	HK2	HK2	HK2	HK2	HK2	HK2	HK2	HK2	HK2
IVT	IVT	IVT	IVT	IVT	IVT	IVT	IVT	IVT	IVT	IVT	IVT
POS	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS
gD	gD	gD	gD	gD	gD	gD	gD	gD	gD	gD	gD
NRT	NRT	NRT	NRT	NRT	NRT	NRT	NRT	NRT	NRT	NRT	NRT
POS	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS	POS
NTC	NTC	NTC	NTC	NTC	NTC	NTC	NTC	NTC	NTC	NTC	NTC

- HK1, HK2 High and Low Level Housekeeping Genes
- IVT RNA Control spiked-in from C-02
- POS Positive PCR Control ± Template
- gD Detects Genomic DNA Contamination
- NRT No Reverse Transcription Control
- NTC No Template Control



Supporting Existing Customers

- ❖ All new customers: PCR Array System 2.0
- ❖ Existing Customers:
 - ❖ New features more useful than old
 - ❖ New features help with technical support issues
- ❖ If old features preferred,
 - ❖ Cannot compare results with System 2.0
 - ❖ Finish project in 6 months – discontinuation of System = C-01, master mix, PCR Array 1.0 (APM,APM,APR)



Typical Customer Features

- ❖ Interested in real-time PCR gene expression profiling
- ❖ Focusing on a pathway or disease state
- ❖ Limited time resources necessary for primer design and optimization

- ❖ NEW FEATURES:
- ❖ High-volume core facilities needing consistency controls
- ❖ 384-well: Industrial accounts
- ❖ Naïve users needing help with RNA QC



CASE™ Kits


Cellular Activation Signaling ELISA




SuperArray Bioscience Corporation
Heather Fox-Brashears
Technical Support



Where is CASE™ used in Gene Expression Research?



ELISA-based assays for *in cell*
Detection of Protein Activation




What is CASE™?

- ❖ **CASE™** = **C**ellular **A**ctivation of **S**ignaling **E**LISA
- ❖ A cell-based assay designed to quantify the level of phosphorylation of a specific protein
- ❖ Includes complete antibody-based detection system
- ❖ Monitors activation of signal transduction pathways
- ❖ Determine whether pathway or protein is activated or inhibited by experimental conditions



How It Works

All experimental steps take place in a single 96 well cell culture plate.



Simply compare relative amounts of target phosphorylated protein to target total protein among your samples to identify the effects of your experimental treatment.

How Does the CASE™ Kit Work?



Benefits of CASE Kits

❖ Quantitative:

Cell-based ELISA determines level of total and phosphorylated form of a protein at same time

❖ Simple:

Easy, quantitative, and non-radioactive protocol with minimal hands-on time

❖ No Extractions:

Directly measure protein phosphorylation state on cells fixed in a 96-well culture plate



Protein Activation: Old Techniques

- ❖ Western blotting
 - ☹ requires cell lysate preparation, immunoprecipitation, and can lead to loss of activated protein
 - ☹ laborious multi-step, time consuming process (low throughput)
 - ☹ separate blots required for total / phosphorylated forms of protein
- ❖ *In vitro* kinase assays
 - ☹ requires protein extraction
 - ☹ requires lots of material
 - ☹ sometimes requires radio-labeling



Protein Activation: New CASE Kits

❖ CASE Benefits:

- ☺ No protein extraction or lysate preps – Cell-based
(maintains native activated protein levels)
- ☺ 96-well plate format – high-throughput compatible
- ☺ Simple ELISA-based procedure, non-radioactive
- ☺ Simultaneous detection of total & phosphorylated forms of protein
- ☺ Little material required – minimal hands-on time
- ☺ Quantitative results



CASE Kits Currently Available

<u>Protein</u>	<u>Site</u>	<u>Species</u>	<u>Catalog #</u>	<u>Inhibitor</u>
⌘⌘ AKT	S473	human, mouse	FE-001	FA-002
⌘⌘ ATF2	T69/T71	mouse	FE-020	N/A
⌘⌘ BAD	S112	human, mouse	FE-021	N/A
⌘⌘ BCR	Y177	mouse	FE-019	N/A
⌘⌘ EGFR	Y845	human, mouse	FE-013	N/A
⌘⌘ ErbB2	Y877	human, mouse	FE-012	N/A
⌘⌘ ERK1/2	T202/Y204	human, mouse	FE-002	FA-003
⌘⌘ IκBa	S32	human	FE-008	FA-006
⌘⌘ JNK	T183/Y185	human, mouse	FE-004	FA-005
⌘⌘ JUN	S73	human, mouse	FE-009	FA-005



CASE Kits Currently Available

<u>Protein</u>	<u>Site</u>	<u>Species</u>	<u>Catalog #</u>	<u>Inhibitor</u>
⊞ NFκB p65	S276	human	FE-007	FA-006
⊞ NFκB p65	S468	human	FE-006	FA-006
⊞ NFκB p65	S536	human	FE-005	FA-006
⊞ p38	T180/Y182	human, mouse	FE-003	FA-004
⊞ p53	S9	human	FE-014	N/A
⊞ p53	S15	human	FE-015	N/A
⊞ p53	S37	human	FE-014	N/A
⊞ PIK3R1	YxxM	human, mouse	FE-010	N/A
⊞ SRC	Y418	human, mouse	FE-011	N/A
⊞ STAT3	Y705	human, mouse	FE-018	N/A
⊞ STAT3	S727	human, mouse	FE-017	N/A



Inhibitors Currently Available

<u>Catalog #</u>	<u>Product</u>	<u>Inhibits</u>
⚡FA-002	LY294002	PI3K
⚡FA-003	U0126	MEK1, 2
⚡FA-004	SB202190	p38 MAPK
⚡FA-005	SP600125	JNK1, 2, 3
⚡FA-006	Bay11-7085	IKBKA/B

Each stock inhibitor will serve > 200 assays




A photograph of a long, straight asphalt road stretching into the distance. The road has a double yellow line down the center and white lines on the sides. The landscape is arid with dry grass and shrubs. In the background, there are dark, rugged mountains under a sky with scattered white clouds. The overall scene conveys a sense of a long journey or a clear path forward.

Focus on
your Pathway!

SureSilencing™ shRNA



Genome-Wide, Plasmid-Based RNA Interference For Human, Mouse, and Rat



Technical Sales Training



SureSilencing™ shRNA Plasmids



- Genome-Wide Availability
- Still new and emerging field **but ...**
- Potential Success Story
- Sales Tools
- Important Technical Support Tips




SureSilencing™ shRNA Introduction

- ❖ Knock down gene expression to determine gene function
- ❖ Pre-designed shRNA specific for a given target gene
- ❖ Unique and Powerful Combination of:
 - ❖ Promega GeneClip™ U1 Hairpin Cloning System
 - ❖ SuperArray's experimentally verified design algorithm
- ❖ For every **human, mouse, or rat gene**, four (4) pre-designed shRNA are provided on separate plasmids.
- ❖ All genes also include one negative control shRNA, a scrambled artificial sequence with no sequence identity to either genome.



How RNA Interference Works Experimentally in Mammalian Cells



Benefits of SureSilencing™ shRNA

❖ Select or Enrich:

☺ Plasmids available with either Neomycin Resistance or GFP:
Enables selection for stable transfectants and the study of the long-term effects of gene suppression or enrichment of transient transfectants for short term studies

❖ GUARANTEED:


☺ Knock down expression of any targeted by at least 70 percent!

❖ Convenient & Cost-Effective:


☺ Use standard plasmid-based and lipid-mediated transfection methods. Plasmids provide a renewable source of RNA Interference.



Promega Corporation's GeneClip™ U1 Hairpin Cloning System—hMGFP GeneClip™ U1 Hairpin Cloning System—Neomycin



Promega Catalog Number C8780



Promega Catalog Number C8790



SureSilencing™ shRNA Features

❖ Promega GeneClip™ Plasmid Backbone

- ❖ Bacterial origin of replication and ampicillin-resistance marker
- ❖ Choice of one of two mammalian markers is available and included on all plasmids per gene.
 - ❖ Neomycin for stable transfections
 - ❖ GFP for transient transfections
- ❖ U1 promoter, transcribed by RNA Polymerase II
 - ❖ Normally transcribes mRNA
 - ❖ Provides moderate shRNA expression level
 - ❖ Minimizes non-specific off-target and toxic side effects



SureSilencing™ shRNA Features

❖ SuperArray's shRNA Design Algorithm

- ❖ Experimentally verified computer algorithm insures gene-specificity and efficacy.
- ❖ Minimizes off-target effect with an advanced specificity search in addition to BLAST.
- ❖ Includes filters for many of the chemical and sequence properties of shRNA known to be important for activity.
- ❖ Effective (>70%) knock down determined by rigorous real-time RT-PCR assay.
- ❖ Pre-design performed for every human, mouse, and rat gene based on the latest annotation from the NCBI.



TYPICAL CUSTOMER FEATURES

- ❖ **Users tired of repeated siRNA purchases**
- ❖ **Using stable (transformed, 2ndary) but ...**
More-difficult-to-transfect cell lines
 - Not using primary or macrophages or doing animal work
 - Recommend viral delivery system instead
- ❖ **Stable OR Transient transfection applications (50:50)**
- ❖ **Cell biologists with less molecular experience**
 - ❖ Purchase GFP and Validate by Western
 - ❖ Difficult to convince of need for real-time PCR

